


Catalog# BP-50168

Trastuzumab-MMAF

Description:

Trastuzumab-MMAF is an antibody-drug conjugate (ADC) of Anti-Her2 human monoclonal antibody Trastuzumab with MMAF (Monomethylauristatin F) through a non-cleavable linker (MC). MMAF is a synthetic antineoplastic agent and a potent inhibitor of tubulin polymerization which applied as a cytotoxic component of the ADCs. The DAR of Trastuzumab-MMAF is $n = 3.5 - 4.0$. This product is used for research only. It can be used as a reference compound. The Trastuzumab-MMAF has the following chemical structure:



Product Details	
Species Reactivity	Human
Published Species	Human
Host/Isotype	Cell culture/ IgG1
Class	Monoclonal
Type	Antibody-drug conjugate
Conjugate	Trastuzumab conjugated with MC-MMAF
DAR	3.5 – 4.0
Immunogen	Human Her2
Form	Liquid
Concentration	1 mg/ml
Purification	Size exclusive column
Storage buffer	50 mM Sodium Phosphate, 0.05% Tween-20, pH7.4
Storage conditions	4°C for short time, -20°C or -80°C for long time.