


Catalog# BP-50188

Sacituzumab govitecan

Description:

Sacituzumab govitecan (IMMU-132) is an anti-human Trop-2 antibody-drug conjugate (ADC) for targeting delivery of the toxin component SN-38. Sacituzumab govitecan shows anticancer activity. It has a molecular weight of 160KDa. The DAR is $n = 7.6$. This product is for research use only. The ADC has the following chemical structure:



Product Details	
CAS No.	1491917-83-9
Species Reactivity	Human
Type	Antibody-drug conjugate
ADC antibody	Sacituzumab (anti-Trop-2)
Target	Human Trop-2
ADC cytotoxin (drug)	SN38
Purity	98.93%
DAR	7.6
Molecular Weight	160kDa
Formulation	10 mg/mL in PBS buffer pH7.4
Purification	Size exclusive column
Storage conditions	Store the undiluted solution at -20°C in the dark to avoid freeze-thaw cycles.