

Catalog# BP-50201

## Telisotuzumab vedotin

### Description:

Telisotuzumab vedotin (Teliso-V, ABBV-399) is an antibody-drug conjugate (ADC) composed of the anti-c-Met humanized monoclonal antibody coupled to the cytotoxic monomethyl auristatin E (MMAE) through a mc-val-cit-PABC type linker. Telisotuzumab vedotin can be used for research on advanced solid tumours. The DAR is  $n = \sim 3.71$ . This product is for research use only.

| Product Details      |                                                                                |
|----------------------|--------------------------------------------------------------------------------|
| CAS No.              | 1714088-51-3                                                                   |
| Species Reactivity   | Human                                                                          |
| Type                 | Antibody-drug conjugate                                                        |
| ADC antibody         | Telisotuzumab (Anti-cMet)                                                      |
| Target               | c-Met                                                                          |
| ADC cytotoxin (drug) | MMAE                                                                           |
| Purity               | 99.1%                                                                          |
| DAR                  | $\sim 3.71$                                                                    |
| Molecular Weight     | N/A                                                                            |
| Formulation          | 1x PBS                                                                         |
| Purification         | Size exclusive column                                                          |
| Storage conditions   | Store the undiluted solution at -20°C in the dark to avoid freeze-thaw cycles. |