

Catalog# BP-50214

Cantuzumab ravtansine

Description:

Zilovetamab vedotin (VLS-101) is a novel antibody-drug conjugate comprising the humanized monoclonal antibody zilovetamab and the anti-microtubule cytotoxin monomethyl auristatin E (MMAE). Zilovetamab vedotin binding to tumor cell ROR1 results in rapid internalization, trafficking to lysosomes, antibody-agent conjugate cleavage, and monomethyl auristatin E release. Zilovetamab vedotin induces apoptosis. Zilovetamab vedotin can be used in research of cancer. This product is for research use only.

Product Details	
CAS No.	868747-45-9
Species Reactivity	Human
Type	Antibody-drug conjugate
ADC antibody	Cantuzumab (anti-huC242)
Target	Humam huC242
ADC cytotoxin (drug)	DM4
Purity	99%
DAR	N/A
Molecular Weight	N/A
Formulation	1xPBS
Purification	Size exclusive column
Storage conditions	Store the undiluted solution at -20°C in the dark to avoid freeze-thaw cycles.