

Catalog# BP-50223

Anti-DM1/DM4 antibody

Anti-DM-1&DM-4 Antibody is a Mouse monoclonal antibody. This antibody has been shown to work in applications such as: PK, PD, Immunoassay and ELISA. Maytansinoid is the highly toxic payload for antibody drug conjugates. This product is used for research only.

Product Details	
Species Reactivity	Mouse
Source	mammalian cell
Isotype	Mouse IgG
Class	Monoclonal
Conjugate	Unconjugated
Immunogen	Maytansinoid (DM1/DM4)
Purity	≥95% (SDS-PAGE)
Protein Concentration	1 mg/ml
Formulation	Lyophilized from sterile PBS, PH 7.4
Storage conditions	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.
Application	Competitive immunoassay validation (Competitive ELISA) and other Immunoassay, PK & PD assay for DM1/DM4 payload of Antibody-drug Conjugate (ADC).