

Catalog# BP-50224

Anti-DXD antibody

Anti-DXD&Exatecan Antibody is a Mouse monoclonal antibody. This antibody has been shown to work in applications such as: PK, PD, Immunoassay and ELISA. DXd/Exatecan is the highly toxic payload for antibody drug conjugates. This product is used for research only.

Product Details	
Species Reactivity	Mouse
Source	mammalian cell
Isotype	Mouse IgG1
Class	Monoclonal
Conjugate	Unconjugated
Immunogen	Deruxtecan (DXd) / Exatecan
Purity	≥95% (SDS-PAGE)
Protein Concentration	1 mg/ml
Formulation	Lyophilized from sterile PBS, PH 7.4
Storage conditions	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.
Application	Competitive immunoassay validation (Competitive ELISA) and other Immunoassay, PK & PD assay for DXd & Exatecan payload of Antibody-drug Conjugate (ADC).