

Catalog# BP-50225

Anti-SN38 antibody

Anti-SN38 Antibody is a Mouse monoclonal antibody. This antibody has been shown to work in applications such as: PK, PD, Immunoassay and ELISA. SN-38 is the highly toxic payload for antibody drug conjugates. This product is used for research only..

Product Details	
Species Reactivity	Mouse
Source	mammalian cell
Isotype	Mouse IgG1
Class	Monoclonal
Conjugate	Unconjugated
Immunogen	7-Ethyl-10-hydroxycamptothecin (SN-38)
Purity	≥95% (SDS-PAGE)
Protein Concentration	1 mg/ml
Formulation	Lyophilized from sterile PBS, PH 7.4
Storage conditions	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.
Application	Competitive immunoassay validation (Competitive ELISA) and other Immunoassay, PK & PD assay for SN-38 payload of Antibody-drug Conjugate (ADC).