

Catalog# BP- 50517

Tezepelumab Biosimilar (Human TSLP Monoclonal Antibody)

Tezepelumab (anti-TSLP) is human monoclonal antibody (IgG2 λ) that binds specifically to TSLP, blocking it from interacting with its heterodimeric receptor. Tezepelumab can be used for the research of severe, uncontrolled asthma.

Product Details	
CAS No.	1572943-04-4
Species Reactivity	Human
Source	Mammalian cells
Isotype	Human IgG2 lambda
Class	Monoclonal
Type	Antibody
Clone	Tezepelumab Biosimilar
Conjugate	Unconjugated
Immunogen	Human TSLP protein
Purity	>95%
Molecular Weight	144.6 kDa
Protein Concentration	1 mg/ml
Formulation	0.2 μ M filtered PBS solution, pH 7.4
Storage conditions	4°C for short time, -20°C or -80°C for long time.