

Catalog# BP-50527

Risankizumab Biosimilar, Human IL-23A (p19) Monoclonal Antibody

Risankizumab (BI 655066) is a humanized IgG monoclonal antibody, targeting IL-23 p19 subunit ($K_d < 10$ pM) and inhibiting IL-17 production induced by human IL-23 in mouse splenocytes ($IC_{50} = 2$ pM). Risankizumab can be used to research immunological and inflammatory disorders such as psoriasis vulgaris, psoriatic arthritis, generalized pustular psoriasis and erythrodermic psoriasis.

Product Details	
CAS No.	1612838-76-2
Species Reactivity	Human
Source	Mammalian cells
Isotype	Human IgG1 kappa
Class	Monoclonal
Type	Antibody
Clone	Risankizumab Biosimilar
Conjugate	Unconjugated
Immunogen	The human IL-23 alpha protein
Purity	>95%
Molecular Weight	145.92 kDa
Protein Concentration	1 mg/ml
Formulation	0.2 μ M filtered PBS solution, pH 7.4
Storage conditions	4°C for short time, -20°C or -80°C for long time.