

Catalog# BP- 50560

Vandortuzumab Biosimilar (Anti-STEAP1 Monoclonal Antibody)

Vandortuzumab is a CHO-expressed humanized antibody targeting STEAP1. Vandortuzumab has a huIgG1 type heavy chain and a huk type light chain, with a predicted molecular weight (MW) of 148.74 kDa. The isotype control for Vandortuzumab can be referenced as Human IgG1 kappa, Isotype Control.

Product Details	
CAS No.	N/A
Species Reactivity	Human
Source	Mammalian cells
Isotype	Human IgG1 kappa
Class	Monoclonal
Type	Antibody
Clone	Vandortuzumab Biosimilar
Conjugate	Unconjugated
Immunogen	Human STEAP1 protein
Purity	>95%
Molecular Weight	145.5 kDa
Protein Concentration	1 mg/ml
Formulation	0.2 µM filtered PBS solution, pH 7.4
Storage conditions	4°C for short time, -20°C or -80°C for long time.