

Catalog# BP- 50570

Felvizumab Biosimilar (RSV Monoclonal Antibody)

Felvizumab Biosimilar uses the same protein sequences as the therapeutic antibody felvizumab. Felvizumab is a humanized monoclonal antibody against respiratory syncytial virus.

Product Details	
CAS No.	N/A
Species Reactivity	Human
Source	Mammalian cells
Isotype	Human IgG1 kappa
Class	Monoclonal
Type	Antibody
Clone	Felvizumab Biosimilar
Conjugate	Unconjugated
Immunogen	Human Respiratory Syncytial Virus
Purity	>95%
Molecular Weight	~150 kD
Protein Concentration	1 mg/ml
Formulation	0.2 µM filtered solution, pH 7.4
Storage conditions	4°C for short time, -20°C or -80°C for long time.